

DRVISION TECHNOLOGIES

AIVIA 5.0.1

April 28th 2017, Bellevue, WA, USA

CHANGE LOG 5.0.0 ⇒ 5.0.1

3D

- Improved speed of dendrite editing in Neuron Editor
- Added heatmap as a coloring mode
- Improved mesh generation of somas when importing from SWC files
- Fixed width, height, and depth measurement values for Neuron and Soma measurements
 - Relocated neuron-related measurements from Soma to Neuron measurements
- Fixed case in Analysis Cube causing errors to occur
- Improved center line accuracy in generation of meshes

2D

- Fixed display issue that could cause tracks to disappear in certain cases

Charts

- Improved caching of axis selectors to be quicker after the first time they have been generated

General

- Improved support for long names while exporting Excel data
- More icon changes